

NR 739

Date: 26-May-04

Application Data Sheet

Hydrochloric acid, Iron(II) & Iron(III)

Matrix

Etching bath.

Principle

Iron(II) is determined by titration with potassium dichromate in sulfuric acid medium.
 $6 Fe^{2+} + 14 H^+ + (Cr_2O_7)^{2-} ==> 6 Fe^{3+} + 7 H_2O + 2 Cr^{3+}$

Iron(III) is determined by titration with EDTA.
 $Fe^{3+} + (H_2Y)^{2-} ==> FeY^- + 2 H^+$

HCl is determined by titration with sodium carbonate.
 $2 HCl + Na_2CO_3 ==> 2 NaCl + H_2O + CO_2$

Detection method

	Method:	Detector	Ion:	λ:
Fe2+	Titration - Redox	Pt electrode	n.a.	n.a.
Fe3+	Titration - Complexometric	Pt electrode	n.a.	n.a.
HCl	Titration - Acid/base	pH Glass electrode	n.a.	n.a.

Specification

	Range	Standard Dev.	Repeatability	Inaccuracy	Analysis time
Fe2+	10 - 200 g/l	< 1 g/l	+/- < 3 g/l	+/- < 3 g/l	30 minutes
Fe3+	1 - 20 g/l	< 0.1 g/l	+/- < 0.3 g/l	+/- < 0.3 g/l	in total
HCl	15 - 250 g/l	< 1.4 g/l	+/- < 4 g/l	+/- < 4 g/l	

Interferences

Other acids and metals.

Reagents

H2SO4 (2 M) 10 ml per analysis
 K2Cr2O7 (0.03 M)
 EDTA (0.02 M)
 Na2CO3 (0.2 M)

Procedure

- clean the analysis vessel with water
- take 0.5 ml of filtered sample and transfer to the analysis vessel with 10 ml of water
- add H2SO4 solution
- perform inflection point titration with K2Cr2O7
- clean the analysis vessel with water
- take 0.5 ml of sample and transfer to the analysis vessel with 15 ml of water
- perform endpoint titration with H2SO4 to pH 1.0
- perform inflection point titration with EDTA
- perform inflection point titration with Na2CO3
- calculate results

Remarks

Sample filtration is needed to prevent clogging of the sample loop valve.

Possible Analyzer

- 2040
- 2016
- 2018 HD
- 2019 HD
- 2019 Special
- 2003 Alert
- 2004 Alert

Typical Wet Part layout

(Other layouts may be realised in order to meet desired criteria, e.g measuring range.)

